

Multiplelabs
Turnkey Laboratories Solutions

FRC-9641S

FRACTION COLLECTOR SYSTEM




- **INTEGRATED COLLECT/DRAIN VALVE**
- **EXCHANGABLE COLLECTION RACKS**
- **FREE POSITION SELECTION**
- **SIMPLE I/O INTERFACE**


►► **ORDER INFORMATION**

» Fraction Collector System / Collector Racks FRC-9641S


Description		Specifications
FRC-3641S	incl. Collector Rack 1	81x 14 ml tubes (diameter 16mm)
	incl. Collector Rack 2	45x 50 ml tubes (diameter 30mm)
	incl. Collector Rack 3	9x 250 ml bottles (diameter 71mm)



Rack 1 - 81x 14ml tubes



Rack 2 - 45x 50ml tubes



Rack 3 - 9x 250ml bottles

►► **Instrument Specifications**

Description	Notes
RS-232	Collect/Drain, Init Position, Next Position, Previous Position, Goto Position, Error Reporting
Digital I/O	Collect/Drain, Init Position, Next Position, Previous Position, Error Output
Wetted Materials	PEEK, PVDF
Max. Flowrate	100 ml/min