



Multiplelabs
Turnkey Laboratories Solutions

FC-9100

Fraction Collector



►► Technical Specifications

Display	Blue LCD display in English
Tubes	120(each max xapacity:12ml and 20ml with a waste tank)
Timing range	1 s-999 h59m 59 s (arbitrary)
Drop range	1 drop -9999 drop (arbitrary)
Pack rang	1-200mv(arbitrary)
Programmable rang	1-10sect(arbitrary)
Moves	Hanging X-Y Tube instead with 12V DC solenoid valve
Operation	English menu - style
Temperature	4-35 °C
Procedures available	1-10modes
Data retention with power outages	Ten years
Tube rack	Rectangular design with a waste tank; free combination (SP-169A, 36A ,16A are optional and scalable available for 36L)