


Multiplelabs
Turnkey Laboratories Solutions

PC-9002

Photocolormeter



►► Technical Specifications

Set Zero (Auto Zero)	By Press Push Switch with Audio response
Wavelength Range	405 - 700 nm
Filters	Built-in 8 Digital Filters >30,000 Hours Life
Filters Wavelength	405nm, 450nm, 480nm, 520nm, 540nm, 578nm, 620nm, 700nm
Display Parameters	O.D., %T, Wavelength
Battery Indicator	Graphical Indication
Detector	Silicon Photo Diode
Accuracy	+ 0.01 O.D.
Stability	± 0.02 O.D. Per Hour
Light Source	Highly Accurate LED Source
Volume	1 ml Minimum
Dimension	275 x 187 x 94mm (Approx.)
Weight	1 Kgs (Approx.)
Power	230 V + 10% AC, 50 Hz
Battery Backup	7.2 V Li-ion Re-Chargeable Battery
Accessories	Set of Test Tube- 5 Nos. Test Tube Holder 1 No. Operating Manual

►► Salient Features

- Digital Wavelength Selection
- Auto Zero Blanking
- Dual Parameter Display